

Runbo

Product name: KN95 Respirator Face Mask
(non-medical, non-surgical)

Model: AKF2002

Standard: GB2626-2006

User Instructions:


1. Failure to follow all instructions and limitations could seriously reduce the effectiveness of this particle filtering half mask and could lead to illness, injury, or death.
2. A properly selected respirator is essential, before occupational use, the wearer must be trained by the employer in the correct use of the respirator in accordance with applicable safety and health standards.
3. This particle filtering half mask does not supply oxygen. Use it in adequately ventilated and where containing sufficient oxygen to support life.
4. Discard this respirator and replace with a new one if:
 - a. Excessive clogging of the particle filtering half mask cause breathing difficulty.
 - b. The particle filtering half mask becomes damaged
5. Leave the contaminated area if dizziness, irritation, or other distress occurs.

Use Limitation:

1. Do not use the respirator or enter or stay in a contaminated area under the following circumstances:
 - a. Atmosphere contains less than 19.5% oxygen.
 - b. If you smell or taste contaminant
 - c. For protection against gases and vapors
 - d. Contaminants or their concentrations are unknown or immediately dangerous to life or health.
 - e. For sandblasting, paint-spray operations and asbestos.
 - f. In explosive atmospheres.
2. Do not modify or misuse the particle filtering half mask
3. Do not use the particle filtering half mask with facial hair or any other conditions that may prevent a good face-seal, the requirements for leakage will be achieved.
4. Particle filtering half mask need to be inspected prior to each use to assure there are no holes in the breathing zone other than punctures around the staples and no damage has occurred. Enlarged holes resulting from ripped or torn filter material around staple punctures are considered damage.
5. This respirator helps protect against certain particulate contaminations but does not eliminate exposure to out the risk of contracting disease or infection. Misuse may result in sickness or death.
6. This particle filtering half mask marked "NR", shall not be used for more than one shift.
7. This particle filtering half mask is not suitable for users who have an ear that has been damaged.

Fitting Instruction:

1. Place finger inside the particle filtering half mask. Bend the nose clip around fingers to form nose shape.
2. Hold the particle filtering half mask in position over the nose and mouth, pull the head harness behind the ears.
3. Use both hands, form the nose clip around the nose and shift the head harness to ensure a positive face seal.
4. Seal check.
 - a. To test the fit of the particle half mask without an exhalation valve, cup both hands over particle filtering half mask and inhale sharply.
 - b. If air flow is felt in the nose area, re-adjust and tighten the nose clip.
 - c. If air flow is felt around the edge of the respirator, adjust the particle half particle filtering half mask/head harness to achieve a better fit.
5. Change the particle filtering half mask immediately if breathing becomes difficult or particle filtering half mask becomes damaged or distorted.
6. Change the particle filtering half mask if proper face seal cannot be achieved.
7. Carefully observe these instructions in an important step in safe respirator use.



Temperature range of conditions



Maximum relative humidity of storage conditions



See information supplied by the manufacturer

Manufacturer: Anshun Health And Medical Technology Co., Ltd.

Address: E6 building industry science and technology park, xixiu district, anshun city, Guizhou province, China.